

Product Name: Human Proinsulin,

Catalog #: AA-0015

Country of Origin: United States of America

Sequence:

FVNQHLCGSHLVEALYLVCGERGFFYTPKTRREAEDLQVGQVELGGGPGAGSLQPLALEGSL
QKRGIVEQCCTSICSLYQLENYCN



UniProt ID: P01308

Counter Ion: Ammonium Bicarbonate

Source: Recombinant peptide expressed in E. Coli

Storage: Lyophilized peptide is shipped at ambient temperature.

**Characterization:
PAGE and HPLC:**



Research Use Only – Not for Use in Humans